LANAFORM®

Relax Mass

Massager cushion

**In-depth shiatsu
massage**

**Massage shiatsu
en profondeur**

**Diepgaande
shiatsu massage**

INTRODUCTION

Thank you for choosing Lanaform's "RELAX MASS", a high quality massage appliance.

"RELAX MASS" is a massage seat cover equipped with two massaging heads to help reduce tension and muscle pain, and soothe tired muscles. Comfortable and improved, the "RELAX MASS" gives a very pleasant 'shiatsu' massage. This appliance is also equipped with a vibration function for the seat and a 'demo' function offering a pleasant blend of the different massage options. In addition, it is operated with a user-friendly control which fits comfortably into the palm of your hand. "RELAX MASS" can be used in the office, at home, ...

PLEASE READ ALL THE INSTRUCTIONS PRIOR TO USING YOUR "RELAX MASS", IN PARTICULAR THESE FEW BASIC SAFETY INSTRUCTIONS:

- Only use this appliance for the purposes described in this manual.
- This appliance is not intended for use by persons (including children) whose physical, sensory or mental capabilities are impaired, or by persons lacking experience or knowledge, unless they have had the benefit of prior supervision and instruction in using the device from a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
- Do not use accessories that are not recommended by LANAFORM® or that are not supplied with this appliance.
- If the power supply lead is damaged, it must be replaced by a similar lead available from the supplier or its after-sales service.
- Do not use this appliance if the plug is damaged, if it is not working correctly, if it has been dropped or if it is damaged or if it has fallen into water. If this is the case, have the appliance inspected and repaired by the supplier or its after-sales department.
- Do not carry this appliance by its supply lead or use this lead as a handle.
- Never use drawing pins or other metal attachments with this appliance.
- Always disconnect the appliance from the mains after use or before cleaning it.
- Should the appliance fall into water, unplug it immediately before retrieving it.
- An electrical appliance should never be left plugged in and unsupervised. Disconnect it from the mains if you are not using it.
- Keep the power lead away from hot surfaces.
- Do not use this product before going to bed. Massage has a stimulating effect which may inhibit sleep.
- Never use this massage appliance while you are asleep.
- Never let anything fall into and never insert anything into one of the openings.
- Never use this appliance in a room where aerosol products (sprays) are used or in a room where oxygen is administered.
- Do not use this appliance under a cover or under a cushion. Excessive heat may cause a fire, electrocution of the individual or other injuries.
- When unplugging the appliance, please make sure that the "On/Off" button is switched off before removing the plug.
- If you are worried about your health, consult a doctor before using this appliance.
- If you feel any pain while using this appliance, stop using it immediately and consult your doctor.
- Do not use this appliance for more than 15 minutes at a time.
- Never use this appliance on a part of the body which is swollen or inflamed or if you have any form of skin rash.
- This product is a non-professional massage appliance designed to soothe tired muscles. Do not use this appliance as a replacement for medical care, it is intended exclusively for family use.

HOW TO USE

1) Details:


- 2 massage heads
- 'Shiatsu' massage
- 3 massage areas (upper back / lower back / full back).
- 1 demo 'function'
- Vibrating seat

2) Description of the "RELAX MASS":

The strap system:

"RELAX MASS" has a system of straps enabling you to fix it to any chair. Simply pass the straps around the chair back, bring them together and fix them to each other using the two adhesive ends.

Control device:

3) Using the "RELAX MASS":

"RELAX MASS" operates via an adapter. First connect the end of the power supply cable to the female input connector located on the right-hand side of the seat. Then plug the adaptor into a power socket. Once connected to the electrical socket, the "RELAX MASS" is ready to use.

- On/Off

Starting or stopping the appliance is controlled by simply pressing the red button (diagram, point 1).

Back:

• 'Demo' mode (diagram, point 2)

This mode consists of a sequence of different massages. This function is valuable since it offers a varied massage.

• Manual mode

You also have the luxury of selecting a specific area or the whole of your back; depending on the desired massage area:

- The whole of your back (diagram, point 3).
- Your lower back (diagram, point 4).
- Your upper back (diagram, point 5).

Seat:

You also have the benefit of a vibrating seat that you can choose to enable or not (diagram, point 6).

Storage

- Place the "RELAX MASS" in its box or in a dry place.
- Avoid it coming into contact with sharp edges and other pointed objects which might cut or pierce the fabric surface.
- In order to avoid any breakage of the power lead, do not wind it around the cushion.
- Do not hang the cushion up by the manual control lead or the power lead.

Cleaning

- Disconnect the cushion from the power socket and wait for a few minutes before cleaning it.
- Clean it with a soft, slightly damp sponge.
- Never let water or any other liquid come into contact with the cushion.
- Do not immerse it in liquid to clean it.
- Never use abrasive detergents, brushes or a thinning agent to clean it.

ADVICE REGARDING THE DISPOSAL OF WASTE

- All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.
- When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law.

LIMITED WARRANTY

LANAFORM® guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.

The LANAFORM® guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM® product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM®'s control.

LANAFORM® may not be held liable for any type of circumstantial, indirect or specific damage.

All implicit guarantees relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.

On receipt, LANAFORM® will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM®'s Service Centre®. Any attempt to maintain this product by a person other than LANAFORM®'s Service Centre will render this guarantee void.

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